manual

Foot massager

Item: 47104

EN

NO

SV

www.rubicson.com
Box 50435 Malmö
Sweden
2021-04-28

Rubicson®

Overview

1. Foot cover
2. Massage heads
3. Massage cushion
4. Heat button
5. Power button
6. Power input
7. Zipper

Specifications

Timer:	15 min
Fabric:	Micro plush, wool
Size:	31x29.5x17 cm
Weight:	1.3 kg

Use

1. Connect the power adapter to the power input (6) and a wall outlet.
2. Choose to use as either foot massage or back massage:

Foot massage

Zip the foot cover (1) onto the massage cushion (3). Sit down in a comfortable chair or couch and place your feet into the foot cover. Do not stand up during the massage.

Back massage

Unzip the foot cover (1). Sit down in a chair or couch with a backrest and place the massager on the part of your back that needs the massage.

3. Press the Power button (5) to turn on the massager.
4. Press the Heat button (4) to turn on heat during the massage. Press again to turn off the heat.

The massager turns off automatically after 15 minutes. Press the Power button to turn it off sooner.

Disconnect the power adapter after use.

Cleaning and maintenance

Clean the massager with a damp cloth if it has become dirty. Clean the rest of the massager with a dry cloth. Do not use any abrasive detergents. These can damage the surface.

Store the massager in a dry environment and at temperatures between 0 and 40 °C.

Safety information

This appliance is not intended for use by persons (including children) with reduced physical, sensory or mental capabilities, or lack of experience and knowledge unless they have been given supervision or instruction concerning use of the appliance by a person responsible for their safety.

Children should be supervised to ensure that they do not play with the appliance.

Use with the provided power adapter only.

If you are unsure whether to use the massager due to a medical condition, consult a doctor.

Do not use the massager for more than 30 continuous minutes. Do not massage one muscle group for longer than 15 minutes.

Do not use the massager if any part of the massager or the power adapter is damaged.

For domestic use only.

Do not immerse the massager in water.

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